

Oxivir® Excel Wipes

1 WIPE - 1 MINUTE



How to Use Oxivir® Excel Wipes



Remove one wipe from pack. Close Lid.



Work from clean to dirty in an S shape. Allow the surface to air dry



Dispose used wipes into a suitable waste bin

READY TO USE PRE-WETTED CLEANING AND DISINFECTING WIPF

Hospital Grade Disinfectant

ARTG No: AUST L 343051 - Oxivir Excel Wipes

Active: Hydrogen Peroxide 0.36%w/w

EFFICACY DATA:

When used as directed, Oxivir Excel Wipes are effective against a wide variety of pathogenic micro-organisms including Viruses, Bacteria and Yeast.

ENVIRONMENTALLY RESPONSIBLE:

Active ingredient breaks down to oxygen and water.



PRODUCT

Oxivir Excel Wipes 24 x 100 wipes (27 x 20cm wipe size)

CODE 70024385

| Bacteria and Yeast | Contact Time |
|----------------------|--------------|
| MRSA | 1 min |
| CRE | 1 min |
| VRE | 1 min |
| Staphylococcus aurei | us 30 sec* |
| Enterococcus hirae | 30 sec* |
| Pseudomonas aerugi | nosa30 sec* |
| Candida albicans | 1 min* |

| Viruses | Contact Time |
|-----------------------|--------------|
| COVID-19 (SARS-CoV-2) | 15 sec |
| Rotavirus | 30 sec |
| Rhinovirus | 30 sec |
| Influenza A H1N1 | 30 sec |
| Herpes Simplex Type 1 | 1 min# |
| Adenovirus Type 5 | 30 sec** |
| Poliovirus Type 1 | 30 sec** |
| Norovirus | 30 sec** |

** Test Method (dirty conditions) EN14476 # Virucidal Carrier Test *Test Method (dirty conditions) EN13727 Yeast (Candida albicans) Test Method (dirty conditions) EN16615



Contact our friendly Customer Service team today
AUSTRALIA - 1800 647 779 - diversey.com.au
NEW ZEALAND - 0800 803 615 - diversey.co.nz